


Product Datasheet

Desmin antibody (orb20095)


Description	Goat polyclonal to Desmin
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Desmin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx. 55kDa band observed in Mouse Skeletal Muscle lysates (calculated MW of 53.5kDa according to NP_001918.3). Recommended concentration: 0.1-0.3 µg/ml.IHC: Paraffin embedded Human Heart. Recommended concentration: 8µg/ml.
Protein Sequence	RDGEVVSEATQQQHE
Clonality	Polyclonal
MW	53.5
Entrez	1674
Dilution Range	ELISA: 1:2000, WB: 0.3-1 µg/ml, IHC-P: 8µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Skeletal ...



Immunohistochemical staining of Human He...



Immunohistochemical staining of Human He...