

Product Datasheet

ALOX15 antibody (orb20099)

Description	Goat polyclonal antibody to ALOX15
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ALOX15
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx. 80kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 74.8kDa according to NP_001131.3). Recommended concentration: 0.2-0.5 µg/ml.
Protein Sequence	HYKTDVAVKDDPE
Clonality	Polyclonal
MW	74.8
Entrez	246
Dilution Range	ELISA: 1:32000, WB: 0.2-0.5 µg/ml, IHC-P: 2 µg/ml
Expiration Date	12 months from date of receipt.



orb20099
(0.2ug/ml)
staining of
nuclear ...