

Product Datasheet

MBNL1 antibody (orb20101)

Description	Goat polyclonal antibody to MBNL1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MBNL1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx. 45kDa band observed in lysates of cell line HeLa (calculated MW of 41.8kDa according to NP_997176.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	SAEHLTSHKYVTQ
Clonality	Polyclonal
MW	41.0; 39.8; 41.8; 33.1; 34.2; 37.0; 37.1
Entrez	4154
Dilution Range	ELISA: 1:4000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20101
(1µg/ml)
staining of
HeLa lysat...