

Product Datasheet

EVL antibody (orb20102)

Description	Goat polyclonal to EVL
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	EVL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx. 50kDa band observed in lysates of cell line Jurkat (calculated MW of 44.8kDa according to NP_057421.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	SQPPNSSEAGRKP
Clonality	Polyclonal
MW	44.8
Entrez	51466
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20102
(1µg/ml)
staining of
Jurkat lys...