


Product Datasheet

MEIS1 antibody (orb20103)

Description	Goat polyclonal to MEIS1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MEIS1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx. 40kDa band observed in Human Brain (Cerebral Cortex) lysates (calculated MW of 43.0kDa according to NP_002389.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	SEDITRSANLTDQ
Clonality	Polyclonal
MW	43
Entrez	4211
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Cerebral ...