


Product Datasheet

GABRB3 antibody (orb20106)

Description	Goat polyclonal antibody to GABRB3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	GABRB3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.IHC: In paraffin embedded Human Cortex shows patchy staining near or at the neuronal cell membrane.. Recommended concentration, 2-4ug/ml.WB: Preliminary experiments in Human Brain (Cerebral Cortex, Frontal Cortex and Cerebellum) lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	EKTAKAKNDRSK
Clonality	Polyclonal
MW	54.1; 54.4
Entrez	2562
Dilution Range	ELISA: 1:8000, WB: 1 µg/ml, IHC-P: 2-4µg/ml
Expiration Date	12 months from date of receipt.



orb20106
(2.5ug/ml)
staining of
paraffin...