

Product Datasheet

ZNF228 antibody (orb20112)

Description

Goat polyclonal to ZNF228

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Target

ZNF228 / ZFP112

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

This antibody is expected to recognize both reported isoforms (NP_001076804.1; NP_037512.3)ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Preliminary experiments gave an approx. 48kDa band in lysates of cell line HepG2 after 2 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 105kDa according to NP_037512.3. The 48kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

EREDDIENSHLKSQ

Clonality

Polyclonal

MW

106; 105

Entrez

7771

Dilution Range

ELISA: 1:4000, WB: 2 µg/ml

Expiration Date

12 months from date of receipt.