



Product Datasheet

EPM2A antibody (orb20119)

Description	Goat polyclonal antibody to EPM2A
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Laforin (isoform a)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx. 38kDa band observed in Human Brain (Cerebellum), Mouse Brain and Rat Brain lysates (calculated MW of 37.2kDa according to NP_005661.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	EATGHTNEMKHTTD
Clonality	Polyclonal
MW	37.2; 35.5
Entrez	7957
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Cerebellu...



Immunohistochemical staining of Human Ce...