

Product Datasheet

VGLL2 Peptide (orb2014221)

Description	VGLL2 Peptide
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
MW	33kDa
Uniprot ID	Q8N8G2
NCBI	NP_872586
Expiration Date	6 months from date of receipt.