

Product Datasheet

SLC26A5 antibody (orb20144)

Description	Goat polyclonal antibody to SLC26A5
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	SLC26A5 / prestin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	ERLHTKDKVPDSIAD
Clonality	Polyclonal
MW	81.3; 74.8; 56.5; 36.9
Entrez	375611
Dilution Range	ELISA: 1:32000
Expiration Date	12 months from date of receipt.