

Product Datasheet

MED15 antibody (orb20146)

Description	Goat polyclonal antibody to MED15
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PCQAP / MED15
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx. 90kDa band observed in Human Lung lysates (calculated MW of 86.8kDa according to NP_001003891.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	PPQSQTQPLVSQAQ
Clonality	Polyclonal
MW	86.8; 82.6
Entrez	51586
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20146
(1µg/ml)
staining of
Human
Lung...