

Product Datasheet

NOA1 antibody (orb20153)

Description

Goat polyclonal antibody to NOA1

Species/Host

Goat

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA

Applications
Target

NOA1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments gave an approx 30kDa band in Human Breast and in Breast cancer lysates after 0.5 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 78.5kDa according to NP_115689.1. The 30kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

ESDYCTAKGSEIDR

Clonality

Polyclonal

MW

78.5

Entrez

84273

Dilution Range

ELISA: 1:32000, WB: 0.5 µg/ml

Expiration Date

12 months from date of receipt.