


Product Datasheet

ARNTL antibody (orb20156)

Description	Goat polyclonal antibody to ARNTL
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	BMAL1 / ARNTL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx. 70+75kDa bands observed in Human Brain (Cerebellum) lysates (calculated MW of 68.7kDa according to NP_001169.3 and 64.1kDa according to NP_001025444.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	REKITTNCYKFKIKD
Clonality	Polyclonal
MW	68.7; 64.1
Entrez	406
Dilution Range	ELISA: 1:8000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Brain (Ce...