


Product Datasheet

NTN1 antibody (orb20164)

Description	Goat polyclonal antibody to NTN1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	Netrin 1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx. 70kDa band observed in Human Brain (Cerebellum), Mouse and Rat Brain lysates (calculated MW of 67.7kDa according to Human NP_004813.1 and 67.8kDa according to Rat NP_446183.1). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	KFQQREKKGKCKKA
Clonality	Polyclonal
MW	67.7
Entrez	9423
Dilution Range	ELISA: 1:32000, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis
of Rat Brain lysat...



Immunohistochemical
staining of Human
In...