

Product Datasheet

ARL4D antibody (orb20179)

Description	Goat polyclonal antibody to ARL4D
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ARL4D
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx. 22kDa band observed in Human Lung lysates (calculated MW of 22.2kDa according to NP_001652.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	KQDQPGALSAAEVEKR
Clonality	Polyclonal
MW	22.2
Entrez	379
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb20179
(0.5ug/ml)
staining of
Human Lu...