

Product Datasheet

UNC5C antibody (orb20184)

Description

Goat polyclonal antibody to UNC5C

Species/Host

Goat

Reactivity

Canine, Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, IF, IHC

Applications
Target

UNC5C

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

IHC: In paraffin embedded Human Cerebral Cortex shows staining of some neuronal bodies plus intense axonal staining. Recommended concentration, 3-5 µg/ml. WB: Preliminary experiments in Human and Mouse Brain lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates

Experiment Notes:

Immunofluorescence: Staining of paraffin embedded endocrine islets (red) in adult Mouse Pancreas. Nuclear staining in blue and ductal cell staining in green. Data obtained from anonymous customer.

Protein Sequence

NCTVSEPTGID

Clonality

Polyclonal

MW


103

Entrez
8633
Dilution Range


ELISA: 1:32000, WB: 1 µg/ml), IHC-P: 3-5 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...