

Product Datasheet

SRX1 antibody (orb20192)

Description	Goat polyclonal to SRX1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SRX1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx. 16kDa band observed in Human Colon and Kidney lysates (calculated MW of 14.3kDa according to NP_542763.1). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	KVQSLVDTIREDPD
Clonality	Polyclonal
MW	14.3
Entrez	140809
Dilution Range	ELISA: 1:32000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb20192
(0.5ug/ml)
staining of
Human Ki...