


Product Datasheet

GRN antibody (orb20199)

Description	Goat polyclonal antibody to GRN
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	Granulin / GRN
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000. WB: Approx. 60kDa band observed in lysates of A549 and HeLa (calculated MW of 63.5kDa according to NP_002078.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	QSKCLSKENATTD
Clonality	Polyclonal
MW	63.5
Entrez	2896
Dilution Range	ELISA: 1:1000, WB: 2 µg/ml, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of paraforma...