


Product Datasheet

E2F7 antibody (orb20204)

Description	Goat polyclonal antibody to E2F7
Species/Host	Goat
Reactivity	Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	E2F7
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 100kDa band observed in Human Brain (Cerebellum), Kidney and Testis lysates (calculated MW of 99.9kDa according to NP_976328.2). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	DLSKQKKFTPERNP
Clonality	Polyclonal
MW	99.9
Entrez	144455
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb20204
(0.1ug/ml)
staining of
Human Ki...