


Product Datasheet

CAV3 antibody (orb20210)

Description	Goat polyclonal antibody to CAV3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CAV3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Skeletal Muscle shows peripheral staining of fibres. Recommended concentration, 3-5µg/ml.WB: Approx. 23kDa band observed in Human, Mouse and Rat Heart lysates (calculated MW of 17.3kDa according to NP_203123.1). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	EHTDLEAQIVKDIH
Clonality	Polyclonal
MW	17.3
Entrez	859
Dilution Range	ELISA: 1:4000, WB: 0.1-0.3 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...



Immunohistochemical analysis of human sk...



Western blot analysis of staining of Pig...