

Product Datasheet

DLX5 antibody (orb20219)

Description	Goat polyclonal antibody to DLX5
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	DLX5
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 30kDa band observed in Human Bone Marrow lysates (calculated MW of 31.5kDa according to NP_005212.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Lee et al, J Biol Chem. 2003 Sep 5;278(36):34387-94; PMID: 12815054). Recommended concentration: 1-3 µg/ml.
Protein Sequence	AYNRVPSATNQPEK
Clonality	Polyclonal
MW	31.5
Entrez	1749
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



orb20219
(1µg/ml)
staining of
Human
Bone...