


Product Datasheet

IFT57 antibody (orb20221)

Description	Goat polyclonal antibody to IFT57
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	HIPPI
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 49kDa band observed in Human and Mouse Brain lysates (calculated MW of 49.1kDa according to Human NP_060480.1 and of 48.8kDa according to Mouse NP_082956.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	EQPQEYDDPNATISN
Clonality	Polyclonal
MW	49.1
Entrez	55081
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb20221 (1ug/ml)
staining of Mouse
Brai...



Immunohistochemical
staining of Human
Te...