


Product Datasheet

TSPAN8 antibody (orb20227)

Description	Goat polyclonal antibody to TSPAN8
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TM4SF3 / TSPAN8
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 24kDa band observed in Human Colon lysates (calculated MW of 26.0kDa according to NP_004607.1). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	QRPCQSYNGKQVYKE
Clonality	Polyclonal
MW	26
Entrez	7103
Dilution Range	ELISA: 1:16000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



orb20227
(0.01ug/ml)
staining of
Human C...