


Product Datasheet

F2R antibody (orb20238)

Description	Goat polyclonal antibody to F2R
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	F2R / PAR1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 48kDa band observed in lysates of cell line HeLa (calculated MW of 47.4kDa according to NP_001983.2). Recommended concentration: 0.5-2 µg/ml.
Protein Sequence	TRARRPESKATNAT
Clonality	Polyclonal
MW	47.4
Entrez	2149
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb20238 (0.5ug/ml)
staining of HeLa lys...



Immunohistochemical
staining of Human
Li...