

Product Datasheet

IL18 antibody (orb20245)

Description	Goat polyclonal antibody to IL18
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	IL18
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 23kDa band observed in Human Placenta lysates and in lysates of cell line HeLa (calculated MW of 22.3kDa according to NP_001553.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	NPPDNIKDTKSDI
Clonality	Polyclonal
MW	22.3; 21.9
Entrez	3606
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20245
(0.1ug/ml)
staining of
HeLa lys...