

Product Datasheet

PON2 Protein (orb2024685)

Description	PON2 Protein
Species/Host	E. coli
Tested Applications	ELISA, IP, WB
Preservatives	PON2 is supplied in PBS and 50% glycerol. Physical appearance: Sterile Filtered clear solution.
Form/Appearance	PON2 is supplied in PBS and 50% glycerol. Physical appearance: Sterile Filtered clear solution.
Storage	Store at 4C if entire vial will be used within 1-2 weeks. Store, frozen at -2°C for longer periods of time. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Application Info: Arylesterase 2 can be used directly as a positive control in Western blotting, ELISA, immunoprecipitation and other immunological experiments. The biological activity of this product has not yet been tested.
Protein Sequence	MGRLLVAVGLLGLIALALLGERLLALRNRLKASREVESVDLPHCHLIKGIEAGSEDID ILPNGLAFFSVGLKFPGLHSFAPDKPGGILMMDLKEEKPRARELRISRGFDLASFN P HGISTFIDNDDTVYLFVVNHPEFKNTVEIFKFEEAENSLHLKTVKHELLPSVNDIT AVGPAHFYATNDHYFSDPFLKYLETYLNLHWANVVYYSPNEVKVVAEGFDSAN GINISPDDKYIYVADILAHEIHVLEKHTNMNLTQLKVELEDTLVDNLSIDPSSGDIW VGCHPNGQKLFVYDPNNPPSSEVLRIQNILSEKPTVTTVYANNGSVLQGSSVASV Y DGKLLIGTLYHRALYCELZ.
Purity	Greater than 95% as determined by SDS-PAGE. Single band on Western Blot.
Uniprot ID	Q15165
NCBI	NP_000296.2
Expiration Date	6 months from date of receipt.