

Product Datasheet

IL17RA antibody (orb20249)

Description	Goat polyclonal antibody to IL17RA
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IL17ra (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 110kDa band observed in Rat Kidney lysates (calculated MW of 96.1kDa according to Mouse NP_032385.1 and of 97.9kDa according to Rat NP_001101353.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	DQEKHGDDSKIN
Clonality	Polyclonal
MW	96.1
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20249
(1ug/ml)
staining of
Rat Kidney...