



Product Datasheet

41888 antibody (orb20255)

Description	Goat polyclonal antibody to 41888
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SEPT6
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 50kDa band observed in Human Testis, Thyroid Gland and MOLT4 cell lysates (calculated MW of 49.7kDa according to NP_055944.2). Recommended concentration: 0.03-1 µg/ml.
Protein Sequence	DEVNAFKQRKTA
Clonality	Polyclonal
MW	48.9; 49.7; 49.2
Entrez	23157
Dilution Range	ELISA: 1:4000, WB: 0.03-1 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Testis ly...



Immunohistochemical staining of Human To...