


Product Datasheet

DCN antibody (orb20265)

Description	Goat polyclonal antibody to DCN
Species/Host	Goat
Reactivity	Bovine, Guinea pig, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Decorin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 38kDa band observed in Human Kidney lysates (calculated MW of 39.7kDa according to NP_001911.1). Recommended concentration: 0.5-2 µg/ml.
Protein Sequence	KISRVDAAASLKGLNN
Clonality	Polyclonal
MW	39.7; 27.4; 23.3;
Entrez	1634
Dilution Range	ELISA: 1:8000, WB: 0.5-2 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb20265
(0.5ug/ml)
staining of
Human Ki...