


Product Datasheet

KEO4 antibody (orb20266)

Description	Goat polyclonal to KEO4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	KEO4 / ERLIN1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Reported variants represent identical protein: NP_006450.2, NP_001094096.1ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 38kDa band observed in Human Placenta lysates (calculated MW of 39.2kDa according to NP_006450.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	EALEPSGENVIQNKE
Clonality	Polyclonal
MW	39.2
Entrez	10613
Dilution Range	ELISA: 1:8000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb20266
(0.5ug/ml)
staining of
Human Pl...