



Product Datasheet

ITGAL antibody (orb20269)


Description	Goat polyclonal antibody to ITGAL
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CD11A / integrin alpha L
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 200kDa band observed in Human Spleen and Tonsil lysates (calculated MW of 129kDa according to NP_002200.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.5-2 µg/ml.
Protein Sequence	YQVRIQPSIHDHN
Clonality	Polyclonal
MW	129; 119
Entrez	3683
Dilution Range	ELISA: 1:16000, WB: 0.5-2 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Spleen Ly...



Immunohistochemical staining of Human Sp...



Immunohistochemical staining of Human Ly...