


Product Datasheet

PRMT5 antibody (orb20273)

Description	Goat polyclonal to PRMT5
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PRMT5
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both isoforms (NP_006100.2; NP_001034708.1).ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 70kDa band observed in Jurkat lysates (calculated MW of 72.7kDa and 71.3kDa according to NP_006100.2; NP_001034708.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	RACREKDRDPEAQ
Clonality	Polyclonal
MW	72.7; 71.3
Entrez	10419
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20273
staining
(2ug/ml) of
Jurkat lys...