


Product Datasheet

CA1 antibody (orb20280)

Description	Goat polyclonal to CA1
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CA1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Reported variants represent identical protein: NP_001729.1, NP_001122303.1, NP_001122302.1, NP_001122301.1.ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 29kDa band observed in Human Liver lysates (calculated MW of 28.9kDa according to NP_001729.1. Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	QAIKTKGKRAP
Clonality	Polyclonal
MW	28.9
Entrez	759
Dilution Range	ELISA: 1:64000, WB: 0.03-0.1 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...