


Product Datasheet

SUMF1 antibody (orb20281)

Description	Goat polyclonal to SUMF1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SUMF1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 37-40kDa band observed in Mouse Pancreas, Mouse Eye and Human Kidney lysates (calculated MW of 40.6kDa according to Human NP_877437.2 and 40.6kDa according to Mouse NP_666049.2). Recommended concentration: 0.5-2 µg/ml.
Protein Sequence	ETLNPKGPPSGKDR
Clonality	Polyclonal
MW	40.6
Entrez	285362
Dilution Range	ELISA: 1:16000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



orb20281
(1ug/ml)
staining of
Mouse Panc...