


Product Datasheet

TREML1 antibody (orb20286)

Description	Goat polyclonal antibody to TREML1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TREML1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 30kDa band observed in Human Bone Marrow lysates (calculated MW of 32.7kDa according to NP_835468.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	CGPAQNPPNNQTPS
Clonality	Polyclonal
MW	32.7
Entrez	340205
Dilution Range	ELISA: 1:2000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb20286
staining
(0.5ug/ml) of
Human Bo...