


Product Datasheet

PDPN antibody (orb20288)

Description	Goat polyclonal to PDPN
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PDPN
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize all isoforms (NP_006465.3; NP_938203.2; NP_001006625.1; NP_001006626.1). ELISA: Peptide ELISA: antibody detection limit dilution 1:1000. WB: Experiments gave bands at approx 25kDa and 37kDa in Human Brain (Cerebellum) lysates after 1 µg/ml antibody staining. The 37kDa band corresponds to earlier findings in literature with different antibodies (Boucherot et al, Nephrol Dial Transplant. 2002 Jun;17(6):978-84; PMID: 12032185). This protein has a calculated MW of 24.9kDa according to NP_006465.3. Recommended concentration: 1-3 µg/ml.
Protein Sequence	EKVDGDTQTTVEKD
Clonality	Polyclonal
MW	24.9; 24.7; 12.4; 12.2
Entrez	10630
Dilution Range	ELISA: 1:1000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20288
(1µg/ml)
staining of
Human Cere...