


Product Datasheet

POLR2G antibody (orb20289)

Description	Goat polyclonal to POLR2G
Species/Host	Goat
Reactivity	Human, Mouse, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	POLR2G
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx 21kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 19.3kDa according to NP_002687.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	EFDPSNPPCYK
Clonality	Polyclonal
MW	19.3
Entrez	5436
Dilution Range	ELISA: 1:1000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Approx
21kDa band
observed in
Human Brai...