

Product Datasheet

Afamin antibody (orb20305)

Description	Goat polyclonal to Afamin
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Afamin / alpha Albumin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 80kDa band observed in Human Brain (Cerebellum and Cerebral Cortex) lysates (calculated MW of 69.1kDa according to NP_001124.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	NKDDRPKDLSLRE
Clonality	Polyclonal
MW	69.1
Entrez	173
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb20305
(0.3ug/ml)
staining of
Human Ce...