


Product Datasheet

WNT9B antibody (orb20309)

Description	Goat polyclonal antibody to WNT9B
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	WNT15 / WNT9B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 38kDa band observed in Human Brain Cerebellum lysates (calculated MW of 39kDa according to NP_003387.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	KRGNKDLRARADA
Clonality	Polyclonal
MW	39
Entrez	7484
Dilution Range	ELISA: 1:8000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb20309
(0.3ug/ml)
staining of
Human Br...



Western blot
analysis of
staining of
Mou...