


Product Datasheet

NUDC antibody (orb20316)

Description	Goat polyclonal antibody to NUDC
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Target	NUDC
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 40kDa band observed in lysates of cell line HeLa (calculated MW of 38.2kDa according to NP_006591.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	PNLGNGADLPNYR
Clonality	Polyclonal
MW	38.2
Entrez	10726
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb20316
(0.3ug/ml)
staining of
HeLa lys...