


Product Datasheet

FERMT3 antibody (orb20318)

Description	Goat polyclonal antibody to FERMT3
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	URP2 / kindlin-3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 75kDa band observed in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 75.4kDa according to NP_113659.3). Recommended concentration: 1-3 µg/ml.
Protein Sequence	ERARGEELDED
Clonality	Polyclonal
MW	75.4
Entrez	83706
Dilution Range	ELISA: 1:2000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20318
(1µg/ml)
staining of
Human Peri...