



Product Datasheet

NFATC2 antibody (orb20329)

Description	Goat polyclonal antibody to NFATC2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	NFATC2 / NFAT1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 110kDa band observed in Human lysates of cell line MOLT4 (calculated MW of 100kDa according to NP_775114.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	QQFRTGLSSPDARYQ
Clonality	Polyclonal
MW	99.8; 100; 97.7
Entrez	4773
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml, IHC-P: 3-6µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of MOLT4 lysate us...



Immunohistochemical staining of Human Sp...