

Product Datasheet

PKD2 antibody (orb20333)

Description	Goat polyclonal to PKD2
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Polycystin 2 / PKD2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 100kDa band observed in lysates of cell line HeLa (calculated MW of 110kDa according to NP_000288.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	ERAKLRREVLGR
Clonality	Polyclonal
MW	110
Entrez	5311
Dilution Range	ELISA: 1:2000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20333
(1ug/ml)
staining of
HeLa lysat...