


Product Datasheet

AAK1 antibody (orb20340)

Description	Goat polyclonal to AAK1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	AAK1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 100kDa band observed in Human Brain (Amygdala and Cerebellum) lysates (calculated MW of 104kDa according to NP_055726.3). Recommended concentration: 1-3 µg/ml.
Protein Sequence	GSPRTSQQNVYNPSE
Clonality	Polyclonal
MW	104
Entrez	22848
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20340
(1µg/ml)
staining of
Human
Amyg...