



Product Datasheet

LIGHT antibody (orb20341)


Description	Goat polyclonal to LIGHT
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Target	LIGHT / CD258
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms (NP_003798.2; NP_742011.1).ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 22kDa band observed in Human Brain (Cerebellum and Frontal Cortex) lysates (calculated MW of 22.4kDa according to NP_742011.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	HEVNPA AHLTGANS
Clonality	Polyclonal
MW	26.4; 22.4
Entrez	8740
Dilution Range	ELISA: 1:64000, WB: 0.03-0.1 µg/ml, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



orb20341 (1ug/ml)
staining of Human
Cere...



WB analysis of
Human Cerebellum
(A) and ...



Immunofluorescence
analysis of HeLa
cell...