

Product Datasheet

NEFM antibody (orb20343)

Description	Goat polyclonal antibody to NEFM
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	NEFM / NF-M
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 170-200kDa band observed in Mouse Brain lysates (calculated MW of 102kDa according to NP_005373.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	QASHITVERKDYLK
Clonality	Polyclonal
MW	102; 59.5
Entrez	4741
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb20343
(0.03ug/ml)
staining of
Mouse B...