


Product Datasheet

PRDM1 antibody (orb20352)

Description	Goat polyclonal antibody to PRDM1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	PRDM1 / MEL1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Tonsil shows textured nuclear staining. Recommended concentration, 3-5µg/ml.WB: Approx 90kDa band observed in lysates of cell lines A431 and A549 (calculated MW of 91.8kDa according to NP_001189.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	DISDNADRLEDVED
Clonality	Polyclonal
MW	91.8; 76.8
Entrez	639
Dilution Range	ELISA: 1:1000, WB: 0.1-2g/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of A431 lysate (35...