

Product Datasheet

HIVEP2 antibody (orb20361)

Description	Goat polyclonal antibody to HIVEP2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	HIVEP2 / MIBP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Preliminary experiments in Human Lymph Node and Human Spleen lysates and in lysates of cell line MOLT4 gave no specific signal but low background (at antibody concentration up to 2 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	TEESKDPSSSEKSQ
Clonality	Polyclonal
MW	269
Entrez	3097
Dilution Range	ELISA: 1:1000, WB: 2 µg/ml
Expiration Date	12 months from date of receipt.