

Product Datasheet

HIP1R antibody (orb20378)

Description	Goat polyclonal to HIP1R
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA
Target	HIP1R/ HIP12
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Preliminary experiments in Human Brain (Amygdala, Hippocampus and Substantia nigra) lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	RQDHQLDKKDG IYP
Clonality	Polyclonal
MW	119
Entrez	9026
Dilution Range	ELISA: 1:32000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.