


Product Datasheet

CEP290 antibody (orb20381)

Description

Goat polyclonal to CEP290


Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, ICC, IHC

Applications

Immunofluorescence
analysis of cell line...

Target

CEP290

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:128000. WB: Preliminary experiments gave an approx 80kDa band in Human Kidney, Testis, and HeLa lysates after 1 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 290kDa according to NP_079390.3. The 80kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

RNSKHLKQQYRAEN

Clonality

Polyclonal

MW

290

Entrez
80184
Dilution Range

ELISA: 1:128000, IHC-P: 5-7 µg/ml

Expiration Date

12 months from date of receipt.